

INT11 antibody



Catalog Number: 111854

Product name

INT11 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

INT11 Rabbit Polyclonal antibody. Positive WB detected in mouse spleen tissue, mouse colon tissue. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 68 kDa

Preparation

This antibody was obtained by immunization of INT11 recombinant protein (Accession Number: NM_017871). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: 1:20-1:200

Validations



mouse spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:111854(INT11 Antibody) at dilution of 1:300



Immunofluorescent analysis of HepG2 cells using Catalog No:111854(INT11 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)