

Integrin alpha-4 antibody



Catalog Number: 111859

Product name

Integrin alpha-4 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Integrin alpha-4 Rabbit Polyclonal antibody.
Positive IHC detected in human tonsillitis tissue.
Positive WB detected in Jurkat cells. Observed molecular weight by Western-blot: 115 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_000885).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol
pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IHC: 1:20-1:200


Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:111859(ITGA4 antibody) at dilution of 1:3000



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:111859(Integrin alpha-4 Antibody) at dilution



of 1:50

Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:111859(Integrin alpha-4 Antibody) at dilution of 1:50