

# Integrin beta-6 antibody



Catalog Number: 111865

## Product name

Integrin beta-6 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

Integrin beta-6 Rabbit Polyclonal antibody.  
Positive IHC detected in human ovary tumor tissue.

## Preparation

This antibody was obtained by immunization of Integrin beta-6 recombinant protein (Accession Number: NM\_000888). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC

## Dilutions


Recommended Dilution:

IHC: 1:20-1:200

## Validations



Immunohistochemical of paraffin-embedded human ovary tumor using Catalog No:111865(ITGB6 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human ovary tumor using Catalog No:111865(ITGB6 antibody) at dilution of 1:50 (under 40x lens)