

INTS10 antibody

Catalog Number: 111870

Product name

INTS10 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

INTS10 Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive WB detected in human kidney tissue, HeLa cells, rat testis tissue. Observed molecular weight by Western-blot: 70-82 kDa

Preparation

This antibody was obtained by immunization of INTS10 recombinant protein (Accession Number: BC001915). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



human kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:111870(INTS10 antibody) at dilution of 1:500



IP Result of anti-INTS10 (IP:Catalog No:111870, 4ug; Detection:Catalog No:111870 1:500) with HeLa cells lysate 1080ug.