

INTS12 antibody



Catalog Number: 111871

Product name

INTS12 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

INTS12 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells. Observed molecular weight by Western-blot: 49-52kd

Preparation

This antibody was obtained by immunization of INTS12 recombinant protein (Accession Number: XM_005263148). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:111871(INTS12 antibody) at dilution of 1:1000