

Jagged 2 antibody

Catalog Number: 111878

Product name

Polyclonal

Jagged 2 antibody

Ig Type

Rabbit IgG

Specificity

Human; other species not tested.

Applications

ELISA, WB

Antibody description

Dilutions

Recommended Dilution:


WB: 1:200-1:2000

Jagged 2 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 150 kDa

Preparation

Validations

This antibody was obtained by immunization of Peptide (Accession Number: NM_002226). Purification method: Antigen affinity purified.



Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:111878(JAG2 Antibody) at dilution of 1:500

Clonality