

# Jagged 2 antibody



Catalog Number: 111878

## Product name

Polyclonal

Jagged 2 antibody

## Ig Type

Rabbit IgG

## Specificity

Human; other species not tested.

## Applications

ELISA, WB

## Antibody description

Jagged 2 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 150 kDa

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_002226). Purification method: Antigen affinity purified.

## Validations


## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:111878(JAG2 Antibody) at dilution of 1:500