

**Product name**

Polyclonal

JIP3 antibody

**Ig Type****Specificity**

Rabbit IgG

Human, Rat, Mouse; other species not tested.

**Applications****Antibody description**

ELISA, WB

JIP3 Rabbit Polyclonal antibody. Positive WB detected in rat brain tissue. Observed molecular weight by Western-blot: 147 kDa


**Dilutions**

Recommended Dilution:

**Preparation**

WB: 1:200-1:2000

This antibody was obtained by immunization of JIP3 recombinant protein (Accession Number: NM\_015133). Purification method: Antigen affinity purified.

**Validations****Formulation**

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

**Storage**

Store at -20°C. DO NOT ALIQUOT

**Clonality**

rat brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111885(JIP3 Antibody) at dilution of 1:600