

# JIP3 antibody



Catalog Number: 111885

## Product name

Polyclonal

JIP3 antibody

## Ig Type

## Specificity

Rabbit IgG

Human, Rat, Mouse; other species not tested.

## Applications

## Antibody description

ELISA, WB

JIP3 Rabbit Polyclonal antibody. Positive WB detected in rat brain tissue. Observed molecular weight by Western-blot: 147 kDa

## Dilutions

Recommended Dilution:


## Preparation

WB: 1:200-1:2000

This antibody was obtained by immunization of JIP3 recombinant protein (Accession Number: NM\_015133). Purification method: Antigen affinity purified.

## Validations

## Formulation



PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

rat brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111885(JIP3 Antibody) at dilution of 1:600

Store at -20°C. DO NOT ALIQUOT

## Clonality