

JPH1 antibody



Catalog Number: 111896

Product name

JPH1 antibody

Immunogen

[Human JPH1 Recombinant protein \(GST tag\)](#)

Specificity

Human; other species not tested.

Antibody description

JPH1 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells. Observed molecular weight by Western-blot: 72 kDa

Preparation

This antibody was obtained by immunization of JPH1 recombinant protein (Accession Number: NM_001317830). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:1000-1:10000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:111896(JPH1 Antibody) at dilution of 1:2000