

Kallistatin antibody



Catalog Number: 111906

Product name

Kallistatin antibody

Specificity

Human; other species not tested.

Antibody description

Kallistatin Rabbit Polyclonal antibody. Positive WB detected in human spleen tissue, human serum tissue. Positive IHC detected in human pancreas cancer tissue. Observed molecular weight by Western-blot: 49 kDa

Preparation

This antibody was obtained by immunization of Kallistatin recombinant protein (Accession Number: NM_006215). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



human spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:111906(SERPINA4 antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:111906(SERPINA4 antibody) at dilution of 1:50 (under 10x lens)