

# KALRN antibody



Catalog Number: 111907

## Product name

KALRN antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

KALRN Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, human kidney tissue, mouse colon tissue. Observed molecular weight by Western-blot: 200kd

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_007064). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111907(KALRN antibody) at dilution of 1:500