

# kazrin antibody



Catalog Number: 111918

## Product name

kazrin antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

kazrin Rabbit Polyclonal antibody. Positive IHC detected in human colon cancer tissue. Positive WB detected in human heart tissue, mouse brain tissue, mouse heart tissue, rat heart tissue. Observed molecular weight by Western-blot: 46 kDa

## Preparation

This antibody was obtained by immunization of kazrin recombinant protein (Accession Number: NM\_001018000). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:111918(RP1-21O18.1 antibody) at dilution of 1:1200



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:111918(RP1-21O18.1 antibody) at dilution of 1:100 (under 10x lens)