

KCNAB1 antibody

Catalog Number: 111924

Product name

KCNAB1 antibody

IHC: 1:20-1:200

IF: 1:20-1:200


Specificity

human, mouse, rat; other species not tested.

Validations

Antibody description

KCNAB1 Rabbit Polyclonal antibody. Positive WB detected in mouse heart tissue. Positive IF detected in Hela. Positive IHC detected in human brain. Observed molecular weight by Western-blot: 68kd




Preparation

This antibody was obtained by immunization of KCNAB1 recombinant protein (Accession Number: NM_172159). Purification method: Antigen affinity purified.

mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:111924(KCNAB1 antibody) at dilution of 1:1500

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal


Immunohistochemical of paraffin-embedded human brain using Catalog No:111924(KCNAB1 antibody) at dilution of 1:100 (under 10x lens)

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, IF



Dilutions

Recommended Dilution:


WB: 1:500-1:5000

Immunohistochemical of paraffin-embedded

KCNAB1 antibody

Catalog Number: 111924

human brain using Catalog No:111924(KCNAB1 antibody) at dilution of 1:100 (under 40x lens)



Immunofluorescent analysis of HeLa cells, using KCNAB1 antibody Catalog No:111924 at 1:50 dilution and Rhodamine-labeled goat anti-rabbit IgG (red). Blue pseudocolor = DAPI (fluorescent DNA dye).