

ITFG1 antibody



Catalog Number: 111945

Product name

Polyclonal

ITFG1 antibody

Ig Type

Specificity

Rabbit IgG

human, mouse, rat; other species not tested.

Applications

Antibody description

ELISA, WB

ITFG1 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue. Observed molecular weight by Western-blot: 68kd

Dilutions

Recommended Dilution:

Preparation

WB: 1:1000-1:10000

This antibody was obtained by immunization of ITFG1 recombinant protein (Accession Number: NM_030790). Purification method: Antigen affinity purified.

Validations

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:111945(ITFG1 antibody) at dilution of 1:2000

