

# KCNJ1 antibody

Catalog Number: 111968

## Product name

KCNJ1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

KCNJ1 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 75-80 kDa

## Preparation

This antibody was obtained by immunization of KCNJ1 recombinant protein (Accession Number: NM\_000220). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:1000

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111968(KCNJ1 antibody) at dilution of 1:2000



IP Result of anti-KCNJ1 (IP:Catalog No:111968, 4ug; Detection:Catalog No:111968 1:300) with mouse brain tissue lysate 3600ug.