

KCNJ3 antibody



Catalog Number: 111970

Product name

KCNJ3 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

KCNJ3 Rabbit Polyclonal antibody. Positive IF detected in HeLa cells. Positive WB detected in mouse brain tissue.

Preparation

This antibody was obtained by immunization of KCNJ3 recombinant protein (Accession Number: NM_002239). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IF, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

Validations



Immunofluorescent analysis of HeLa cells using Catalog No:111970(KCNJ3 Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111970(KCNJ3 Antibody) at dilution of 1:600