

KCNJ5 antibody

Catalog Number: 111971

Product name

KCNJ5 antibody

Specificity

Human; other species not tested.

Antibody description

KCNJ5 Rabbit Polyclonal antibody. Positive WB detected in transfected cells.

Preparation

This antibody was obtained by immunization of KCNJ5 recombinant protein (Accession Number: XM_011542809). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:250 -1:1000

Validations



WB result of KCNJ5 antibody (Catalog No:111971) with transfected cell lysate. His tag antibody (Catalog No:117294) as control.