

# KCNJ5 antibody



Catalog Number: 111971

## Product name

KCNJ5 antibody

## Specificity

Human; other species not tested.

## Antibody description

KCNJ5 Rabbit Polyclonal antibody. Positive WB detected in transfected cells.

## Preparation

This antibody was obtained by immunization of KCNJ5 recombinant protein (Accession Number: XM\_011542809). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:250 -1:1000

## Validations



WB result of KCNJ5 antibody (Catalog No:111971) with transfected cell lysate. His tag antibody (Catalog No:117294) as control.