

KCNK5 antibody



Catalog Number: 111973

Product name

Polyclonal

KCNK5 antibody

Ig Type

Specificity

Rabbit IgG

human, mouse, rat; other species not tested.

Applications

Antibody description

ELISA, WB

KCNK5 Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells. Observed molecular weight by Western-blot: 72kd

Dilutions

Recommended Dilution:

Preparation

WB: 1:200-1:2000

This antibody was obtained by immunization of KCNK5 recombinant protein (Accession Number: NM_003740). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:111973(KCNK5 antibody) at dilution of 1:800