

KCNN3 antibody



Catalog Number: 111976

Product name

KCNN3 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

KCNN3 Rabbit Polyclonal antibody. Positive IF detected in SH-SY5Y cells. Positive WB detected in Transfected HEK-293 cells.

Preparation

This antibody was obtained by immunization of KCNN3 recombinant protein (Accession Number: NM_170782). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IF, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: 1:10-1:100

Validations



Immunofluorescent analysis of SH-SY5Y cells using Catalog No:111976(KCNN3 Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:111976(KCNN3 Antibody) at dilution of 1:350