

# KIAA0513 antibody



Catalog Number: 111983

## Product name

KIAA0513 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

KIAA0513 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, mouse kidney tissue, mouse cerebellum tissue, human testis tissue. Positive IF detected in HepG2. Observed molecular weight by Western-blot: 34kd

## Preparation

This antibody was obtained by immunization of KIAA0513 recombinant protein (Accession Number: NM\_001297766). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IF: 1:10-1:100

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111983(KIAA0513 antibody) at dilution of



Immunofluorescent analysis of HepG2 cells, using KIAA0513 antibody Catalog No:111983 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).