

KIAA0528 antibody

Catalog Number: 111984

Product name

KIAA0528 antibody

IP: 1:200-1:2000

IF: 1:10-1:100

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

KIAA0528 Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive WB detected in HeLa cells. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 115 kDa

Preparation

This antibody was obtained by immunization of KIAA0528 recombinant protein (Accession Number: NM_014802). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IF, IP

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:111984(KIAA0528 antibody) at dilution of 1:1000



Immunofluorescent analysis of HepG2 cells, using KIAA0528 antibody Catalog No:111984 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit



IgG (red).

KIAA0528 antibody



Catalog Number: 111984

IP Result of anti-KIAA0528 (IP:Catalog

No:111984, 4ug; Detection:Catalog No:111984
1:700) with HeLa cells lysate 520ug.