

# KIAA0528 antibody



Catalog Number: 111984

## Product name

KIAA0528 antibody

IP: 1:200-1:2000

IF: 1:10-1:100

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

KIAA0528 Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive WB detected in HeLa cells. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 115 kDa

## Preparation

This antibody was obtained by immunization of KIAA0528 recombinant protein (Accession Number: NM\_014802). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB, IF, IP

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000


## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:111984(KIAA0528 antibody) at dilution of 1:1000



Immunofluorescent analysis of HepG2 cells, using KIAA0528 antibody Catalog No:111984 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit



# KIAA0528 antibody



Catalog Number: 111984

---

IP Result of anti-KIAA0528 (IP:Catalog

No:111984, 4ug; Detection:Catalog No:111984  
1:700) with HeLa cells lysate 520ug.