

KIAA1033 antibody



Catalog Number: 111988

Product name

KIAA1033 antibody

Ig Type

Rabbit IgG

Specificity

Human; other species not tested.

Applications

ELISA, WB

Antibody description

KIAA1033 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 136 kDa

Dilutions


Recommended Dilution:

WB: 1:200-1:1000

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_015275).
Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:111988(KIAA1033 Antibody) at dilution of 1:300

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal