

# KIAA1191 antibody



Catalog Number: 111989

## Product name

KIAA1191 antibody

## Specificity

human, mouse; other species not tested.

## Antibody description

KIAA1191 Rabbit Polyclonal antibody. Positive WB detected in mouse ovary tissue, mouse uterus tissue. Positive IF detected in Hela. Observed molecular weight by Western-blot: 33kd

## Preparation

This antibody was obtained by immunization of KIAA1191 recombinant protein (Accession Number: BC006316). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

## Validations



mouse ovary tissue were subjected to SDS PAGE followed by western blot with Catalog No:111989(KIAA1191 antibody) at dilution of 1:200



Immunofluorescent analysis of Hela cells, using KIAA1191 antibody Catalog No:111989 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).