

KIAA1530 antibody



Catalog Number: 111994

Product name

KIAA1530 antibody

Specificity

Human; other species not tested.

Antibody description

KIAA1530 Rabbit Polyclonal antibody. Positive WB detected in Raji cells. Observed molecular weight by Western-blot: 81 kDa

Preparation

This antibody was obtained by immunization of KIAA1530 recombinant protein (Accession Number: NM_001317935). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:111994(KIAA1530 Antibody) at dilution of 1:300