

# KIAA1598 antibody



Catalog Number: 111995

## Product name

KIAA1598 antibody

## Specificity

Human, Rat; other species not tested.

## Antibody description

KIAA1598 Rabbit Polyclonal antibody. Positive WB detected in rat brain tissue. Observed molecular weight by Western-blot: 60 kDa

## Preparation

This antibody was obtained by immunization of KIAA1598 recombinant protein (Accession Number: BC022348). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



rat brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111995(KIAA1598 Antibody) at dilution of 1:600