

# KIAA1984 antibody



Catalog Number: 111997

## Product name

KIAA1984 antibody

## Specificity

Human; other species not tested.

## Antibody description

KIAA1984 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HepG2 cells. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 45-55 kDa

## Preparation

This antibody was obtained by immunization of KIAA1984 recombinant protein (Accession Number: BC007542). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:20-1:200

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:111997(KIAA1984 Antibody) at dilution of 1:600



Immunofluorescent analysis of (-20oc Ethanol) fixed HepG2 cells using Catalog No:111997(KIAA1984 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)