

KIF3C antibody

Catalog Number: 112014

Product name

KIF3C antibody

IP: 1:200-1:2000

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

KIF3C Rabbit Polyclonal antibody. Positive IHC detected in human gliomas tissue. Positive IP detected in mouse brain tissue. Positive WB detected in Y79 cells, HeLa cells, mouse brain tissue. Observed molecular weight by Western-blot: 100 kDa

Preparation

This antibody was obtained by immunization of KIF3C recombinant protein (Accession Number: NM_002254). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



Y79 cells were subjected to SDS PAGE followed by western blot with Catalog No:112014(KIF3C antibody) at dilution of 1:800



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:112014(KIF3C antibody) at dilution of 1:100 (under 10x lens)




Immunohistochemical of paraffin-embedded human gliomas using Catalog No:112014(KIF3C

KIF3C antibody

Catalog Number: 112014

antibody) at dilution of 1:100 (under 40x lens)



IP Result of anti-KIF3C (IP:Catalog No:112014, 4ug; Detection:Catalog No:112014 1:500) with mouse brain tissue lysate 3440ug.