

KDM5D antibody



Catalog Number: 112030

Product name

KDM5D antibody

Specificity

Human; other species not tested.

Antibody description

KDM5D Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells.

Preparation

This antibody was obtained by immunization of KDM5D recombinant protein (Accession Number: NM_001146706). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:112030(KDM5D Antibody) at dilution of 1:700