

KIAA0368 antibody



Catalog Number: 112047

Product name

KIAA0368 antibody

Immunogen

[Human KIAA0368 Recombinant protein \(His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

KIAA0368 Rabbit Polyclonal antibody. Positive IF detected in HEK-293 cells. Positive WB detected in HEK-293 cells, A549 cells. Observed molecular weight by Western-blot: 205 kDa

Preparation

This antibody was obtained by immunization of KIAA0368 recombinant protein (Accession Number: BC021127). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:20-1:200

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:112047(KIAA0368 Antibody) at dilution of 1:600



Immunofluorescent analysis of HEK-293 cells using Catalog No:112047(KIAA0368 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)