

KIFC3 antibody



Catalog Number: 112054

Product name

KIFC3 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

KIFC3 Rabbit Polyclonal antibody. Positive WB detected in mouse lung tissue, mouse kidney tissue. Positive IP detected in mouse kidney tissue. Observed molecular weight by Western-blot: 80-85 kDa, 95-100 kDa

Preparation

This antibody was obtained by immunization of KIFC3 recombinant protein (Accession Number: XM_047434084). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:112054(KIFC3 antibody) at dilution of 1:500



IP Result of anti-KIFC3 (IP:Catalog No:112054, 4ug; Detection:Catalog No:112054 1:500) with mouse kidney tissue lysate 7200ug.