

KLF14 antibody

Catalog Number: 112079

Product name

KLF14 antibody

WB: 1:500-1:5000

IHC: 1:20-1:200

Immunogen

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

KLF14 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, mouse liver tissue. Positive IHC detected in human skin tissue, human placenta tissue. Observed molecular weight by Western-blot: 40 kDa

Preparation

This antibody was obtained by immunization of KLF14 recombinant protein (Accession Number: NM_001258248). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:112079(SP6 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human skin using Catalog No:112079(SP6 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human skin using Catalog No:112079(SP6

KLF14 antibody



Catalog Number: 112079

antibody) at dilution of 1:100 (under 40x lens)