

# KLF17 antibody



Catalog Number: 112080

## Product name

Polyclonal

KLF17 antibody

## Ig Type

## Specificity

Rabbit IgG

Human; other species not tested.

## Applications

## Antibody description

ELISA, WB

KLF17 Rabbit Polyclonal antibody. Positive WB detected in MCF-7 cells. Observed molecular weight by Western-blot: 43 kDa

## Dilutions

Recommended Dilution:


## Preparation

WB: 1:200-1:1000

This antibody was obtained by immunization of KLF17 recombinant protein (Accession Number: NM\_173484). Purification method: Antigen affinity purified.

## Validations

## Formulation



PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

MCF-7 cells were subjected to SDS PAGE followed by western blot with Catalog No:112080(KLF17 Antibody) at dilution of 1:300