

# KPNA5 antibody

Catalog Number: 112115

## Product name

KPNA5 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

KPNA5 Rabbit Polyclonal antibody. Positive IP detected in Jurkat cells. Positive WB detected in Jurkat cells. Observed molecular weight by Western-blot: 60 kDa

## Preparation

This antibody was obtained by immunization of KPNA5 recombinant protein (Accession Number: NM\_002269). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IP: 1:1000-1:10000

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:112115(KPNA5 antibody) at dilution of 1:2000



IP Result of anti-KPNA5 (IP:Catalog No:112115, 4ug; Detection:Catalog No:112115 1:2000) with Jurkat cells lysate 4000ug.