

KPNA5 antibody



Catalog Number: 112115

Product name

KPNA5 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

KPNA5 Rabbit Polyclonal antibody. Positive IP detected in Jurkat cells. Positive WB detected in Jurkat cells. Observed molecular weight by Western-blot: 60 kDa

Preparation

This antibody was obtained by immunization of KPNA5 recombinant protein (Accession Number: NM_002269). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IP: 1:1000-1:10000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:112115(KPNA5 antibody) at dilution of 1:2000



IP Result of anti-KPNA5 (IP:Catalog No:112115, 4ug; Detection:Catalog No:112115 1:2000) with Jurkat cells lysate 4000ug.