

**Product name**

KRBP antibody

**Specificity**

Human, Mouse; other species not tested.

**Antibody description**

KRBP Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue. Observed molecular weight by Western-blot: 63kd

**Preparation**

This antibody was obtained by immunization of KRBP recombinant protein (Accession Number: NM\_024106). Purification method: Antigen affinity purified.

**Formulation**

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

**Storage**

Store at -20°C. DO NOT ALIQUOT

**Clonality**

Polyclonal

**Ig Type**

Rabbit IgG


**Applications**

ELISA, WB

**Dilutions**

Recommended Dilution:

WB: 1:200-1:1000

**Validations**

mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:112122(KRBP antibody) at dilution of 1:300