

KRT6B antibody



Catalog Number: 112134

Product name

KRT6B antibody

Recommended Dilution:

WB: 1:500-1:5000

Immunogen

IP: 1:200-1:2000

[Human KRT6B Recombinant protein \(GST tag\)](#)

IHC: 1:20-1:200


Specificity

Human, Mouse, Rat; other species not tested.

Validations

Antibody description

KRT6B Rabbit Polyclonal antibody. Positive IP detected in mouse skin tissue. Positive WB detected in mouse skin tissue, A431 cells. Positive IHC detected in human cervix tissue, human oesophagus tissue, human tonsil tissue. Observed molecular weight by Western-blot: 60 kDa




Preparation

This antibody was obtained by immunization of KRT6B recombinant protein (Accession Number: NM_005555). Purification method: Antigen affinity purified.

mouse skin tissue were subjected to SDS PAGE followed by western blot with Catalog No:112134(KRT6B antibody) at dilution of 1:1000

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Immunohistochemical of paraffin-embedded human cervix using Catalog No:112134(KRT6B antibody) at dilution of 1:100 (under 10x lens)

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC, IP


Dilutions

KRT6B antibody

Catalog Number: 112134



Immunohistochemical of paraffin-embedded human cervix using Catalog No:112134(KRT6B antibody) at dilution of 1:100 (under 40x lens)



IP Result of anti-KRT6B (IP:Catalog No:112134, 4ug; Detection:Catalog No:112134 1:800) with mouse skin tissue lysate 3200ug.