

LANCL1 antibody



Catalog Number: 112139

Product name

LANCL1 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

LANCL1 Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer. Positive WB detected in mouse brain tissue, human heart tissue. Observed molecular weight by Western-blot: 45kd

Preparation

This antibody was obtained by immunization of LANCL1 recombinant protein (Accession Number: NM_001136574). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:112139(LANCL1 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:112139(LANCL1 antibody) at dilution of 1:50 (under 10x lens)