

Langerin antibody



Catalog Number: 112140

Product name

Langerin antibody

Specificity

Human; other species not tested.

Antibody description

Langerin Rabbit Polyclonal antibody. Positive WB detected in human skin tissue. Observed molecular weight by Western-blot: 40kd

Preparation

This antibody was obtained by immunization of Langerin recombinant protein (Accession Number: XM_011532874). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



human skin tissue were subjected to SDS PAGE followed by western blot with Catalog No:112140(CD207 antibody) at dilution of 1:500