

Langerin antibody



Catalog Number: 112140

Product name

Polyclonal

Langerin antibody

Ig Type

Specificity

Rabbit IgG

Human; other species not tested.

Applications

Antibody description

ELISA, WB

Langerin Rabbit Polyclonal antibody. Positive WB detected in human skin tissue. Observed molecular weight by Western-blot: 40kd

Dilutions

Recommended Dilution:

Preparation

WB: 1:200-1:2000

This antibody was obtained by immunization of Langerin recombinant protein (Accession Number: XM_011532874). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



human skin tissue were subjected to SDS PAGE followed by western blot with Catalog No:112140(CD207 antibody) at dilution of 1:500