

# LARP4 antibody



Catalog Number: 112147

## Product name

LARP4 antibody

WB: 1:200-1:2000

IP: 1:200-1:1000

## Specificity


Human, Mouse, Rat; other species not tested.

IHC: 1:20-1:200

## Antibody description

LARP4 Rabbit Polyclonal antibody. Positive WB detected in Y79 cells. Positive IP detected in Y79 cells. Positive IHC detected in human brain tissue, human kidney tissue, human lung tissue, human ovary tissue, human spleen tissue, human testis tissue. Observed molecular weight by Western-blot: 81 kDa

## Validations




## Preparation

This antibody was obtained by immunization of LARP4 recombinant protein (Accession Number: BC022377). Purification method: Antigen affinity purified.

Y79 cells were subjected to SDS PAGE followed by western blot with Catalog No:112147(LARP4 antibody) at dilution of 1:300

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal


Immunohistochemical of paraffin-embedded human brain using Catalog No:112147(LARP4 antibody) at dilution of 1:50 (under 10x lens)

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC, IP



## Dilutions

Recommended Dilution:

# LARP4 antibody

Catalog Number: 112147

---

Immunohistochemical of paraffin-embedded human brain using Catalog No:112147(LARP4 antibody) at dilution of 1:50 (under 40x lens)

IP Result of anti-LARP4 (IP:Catalog No:112147, 4ug; Detection:Catalog No:112147 1:300) with Y79 cells lysate 2000ug.

