

# LARS2 antibody

Catalog Number: 112150

## Product name

LARS2 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

LARS2 Rabbit Polyclonal antibody. Positive IP detected in Jurkat cells. Positive WB detected in HepG2 cells, Jurkat cells, MCF7 cells. Observed molecular weight by Western-blot: 102 kDa

## Preparation

This antibody was obtained by immunization of LARS2 recombinant protein (Accession Number: NM\_015340). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

## Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:112150(LARS2 antibody) at dilution of 1:400



IP Result of anti-LARS2 (IP:Catalog No:112150, 4ug; Detection:Catalog No:112150 1:500) with Jurkat cells lysate 2000ug.