

# LCA5 antibody



Catalog Number: 112167

## Product name

LCA5 antibody

## Immunogen

[Human LCA5 Recombinant protein \(His tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

LCA5 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 80 kDa

## Preparation

This antibody was obtained by immunization of LCA5 recombinant protein (Accession Number: NM\_001122769). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:1000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:112167(LCA5 Antibody) at dilution of 1:300