

LCE1B antibody

Catalog Number: 112170

Product name

LCE1B antibody

IHC: 1:20-1:200

IF: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

LCE1B Rabbit Polyclonal antibody. Positive WB detected in A375 cells, mouse skin tissue, rat skin tissue. Positive IHC detected in human skin tissue, human oesophagus tissue. Positive IF detected in A375 cells. Observed molecular weight by Western-blot: 10-12 kDa

Preparation

This antibody was obtained by immunization of LCE1B recombinant protein (Accession Number: NM_178349). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IF

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



A375 cells were subjected to SDS PAGE followed by western blot with Catalog No:112170(LCE1B antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human skin using Catalog No:112170(LCE1B antibody) at dilution of 1:50 (under 10x lens)




Immunohistochemical of paraffin-embedded human skin using Catalog No:112170(LCE1B antibody) at dilution of 1:50 (under 10x lens)

LCE1B antibody

Catalog Number: 112170

antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of A375 cells using Catalog No:112170(LCE1B Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)