

LDHAL6A antibody



Catalog Number: 112182

Product name

LDHAL6A antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

LDHAL6A Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, human brain tissue, Jurkat cells, mouse brain tissue, mouse testis tissue. Positive IF detected in Hela cells. Observed molecular weight by Western-blot: 30-36 kDa

Preparation

This antibody was obtained by immunization of LDHAL6A recombinant protein (Accession Number: BC014340). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:112182(LDHAL6A antibody) at dilution of



Immunofluorescent analysis of Hela cells, using LDHAL6A antibody Catalog No:112182 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).