

LEF1 antibody

Catalog Number: 112191

Product name

LEF1 antibody

IP: 1:200-1:1000

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

LEF1 Rabbit Polyclonal antibody. Positive IHC detected in human thymus tissue, human spleen tissue. Positive IP detected in SW 1990 cells. Positive WB detected in COLO 320 cells, SW 1990 cells. Observed molecular weight by Western-blot: 40-45 kDa

Preparation

This antibody was obtained by immunization of LEF1 recombinant protein (Accession Number: NM_016269). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC, WB, IP

Dilutions


Recommended Dilution:

WB: 1:200-1:1000


Validations



Immunohistochemical of paraffin-embedded human thymus using Catalog No:112191(LEF1 antibody) at dilution of 1:100 (under 40x lens)



Immunohistochemical of paraffin-embedded human spleen using Catalog No:112191(LEF1 antibody) at dilution of 1:100 (under 40x lens)




COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog

LEF1 antibody

Catalog Number: 112191

No:112191(LEF1 Antibody) at dilution of 1:300



IP Result of anti-LEF1 (IP:Catalog No:112191, 4ug; Detection:Catalog No:112191 1:300) with SW 1990 cells lysate 2400ug.