

LENG1 antibody



Catalog Number: 112193

Product name

LENG1 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

LENG1 Rabbit Polyclonal antibody. Positive IHC detected in human ovary tumor. Positive WB detected in human skin tissue. Observed molecular weight by Western-blot: 31kd

Preparation

This antibody was obtained by immunization of LENG1 recombinant protein (Accession Number: NM_024316). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



human skin tissue were subjected to SDS PAGE followed by western blot with Catalog No:112193(LENG1 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human ovary tumor using Catalog No:112193(LENG1 antibody) at dilution of 1:50 (under 10x lens)