

LILRA4 antibody



Catalog Number: 112224

Product name

LILRA4 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

LILRA4 Rabbit Polyclonal antibody. Positive WB detected in L02 cells, HL-60 cells, mouse spleen tissue. Observed molecular weight by Western-blot: 48kd

Preparation

This antibody was obtained by immunization of LILRA4 recombinant protein (Accession Number: BC109198). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



L02 cells were subjected to SDS PAGE followed by western blot with Catalog No:112224(LILRA4 antibody) at dilution of 1:500